



PIPETTE HOLDERS

Pressure Systems for Ejection of
Picoliter Volumes in Cell Research



P/N 050-00XX-110-1*
Standard straight holder with a
0.080 inch diameter by 1.5 inch
long mounting rod.



P/N 050-00XX-130-1*
Standard straight holder with a
.186 inch diameter by 3.5 inch
long mounting rod.

*XX is replaced by the capillary glass' OD. We have 1.0, 1.2, 1.5,
and 2.0mm options available (Example: P/N 050-0010-110-1 is for a
1.0mm OD capillary glass).